



Accurate & Speedy  
Analyzer  
for Professional

# Labona Check™ **A1c** HbA1c Analyzer




# Labona Check™ A1c

## HbA1c Analyzer



### »» Performance



### »» Features & Benefits

- Auto guide operation
- Large LCD touch-screen
- Automatic tray insertion
- External power source
- Internal printer
- Simple operation step
- User-friendly technology


### »» Specifications

Methodology	Boronate affinity
Assay Time	3 min
Measurement Time	7 sec.
Measuring Method	Optical Reflectance
Temperature Range	10–40 °C
Sample Volume	Whole Blood, 5µl
Measurement Range	3–15%

### »» Test Kit



### »» Accessories (Sold separately)



### »» How to use



Blood sample preparation



Hemoglobin precipitation



Apply the sample and washing solution



Initiate automatic tray



Insert the cartridge



Read the result